PRINCIPIO

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Definition, Diagnosis and Classification of Diabetes Mellitus and its Complications. World Health Organization (1999). Screening for type 2 Diabetes. Clinical Diabetes vol.18 NO. 2 Spring (2000). Tietz, NW., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012). Young, D.S.; Effects of drugs on Clinical Laboratory Tests. 4th Ed. AACCC Press. 2000.

PRINCIPLE

The action mechanism of glucose solutions for the glucose metabolism diagnosis is based on the normal physiological and biochemical processes of the organism. This glucose solution intake cause a glucose overload in the blood from which the physiological capability of the organisms to regulate this glycaemia is evaluated.

The Dextrose Solution GLYCOSULL has been carefully designed to offer a safe and a practical product to standardize glucose metabolism tests. They have been formulated to be pleasant for ingestion and to avoid the frequent nausea caused by these types of preparations. They are products with natural citrus flavourings and do not contain artificial colourings. The additives contained in the product do not have any contraindication according to the current food regulations in the European Union. They have been manufactured following strict physicochemical and microbiological quality controls, under international standards, so batch-to-batch consistency can be ensured.

They are presented in three different concentrations, 50, 75 and 100 g in order to adjust the exact glucose intake according to diagnosis to be carried out and the patient characteristics. Usually the 50g is used for screening gestational diabetes; the one of 75g for the diagnosis of Diabetes mellitus and the one of 100g for the diagnosis of gestational diabetes.

DIAGNOSTIC USE

Glucose solution for oral administration, intended for the performance of diagnostic tests related to glucose metabolism, for glucose tolerance test in the evaluation of diabetes mellitus and related diseases.

Single test result can not be used to make a clinical diagnosis. It should integrate clinical and laboratory data. It must be administered by a healthcare professional and under proper supervision.

REAGENTS

GLYCOSULL ORANGE 50 g	30x200mL	Ref. 99 76 85
GLYCOSULL ORANGE 75 g	30x200mL	Ref. 99 47 32
GLYCOSULL ORANGE 100 g	30x200mL	Ref. 99 61 36
GLYCOSULL LEMON 50 g	30x200mL	Ref. 99 52 71
GLYCOSULL LEMON 75 g	30x200mL	Ref. 99 94 89
GLYCOSULL LEMON 100 g	30x200mL	Ref. 99 59 03

REAGENT COMPOSITION

GLYCOSULL ORANGE 50 g

Citric acid anhydrous (E-330)	0.25 %
Preservative (E-202)	0.1 %
D- Glucose	25 %
Orange flavour	0.1 %

GLYCOSULL ORANGE 75 g

Citric acid anhydrous (E-330)	0.25 %
Preservative (E-202)	0.1 %
D- Glucose	37.5 %
Orange flavour	0.1 %

GLYCOSULL ORANGE 100 g

Citric acid anhydrous (E-330)	0.25 %
Preservative (E-202)	0.1 %
D- Glucose	50 %
Orange flavour	0.1 %

GLYCOSULL LEMON 50 g

Citric acid anhydrous (E-330)	0.25 %
Preservative (E-202)	0.1 %
D- Glucose	25 %
Lemon flavour	0.05 %

GLYCOSULL LEMON 75 g

Citric acid anhydrous (E-330)	0.25 %
Preservative (E-202)	0.1 %
D- Glucose	37.5 %
Lemon flavour	0.05 %

GLYCOSULL LEMON 100 g

Citric acid anhydrous (E-330)	0.25 %
Preservative (E-202)	0.1 %
D- Glucose	50 %
Lemon flavour	0.05 %

All solutions have been prepared from conventional corn dextrose that does not contain wheat gluten (less than 10ppm). There is no risk of cross contamination with gluten residues or products that could contain them in the manufacturing process. Glycosull may contain traces of fructose (<0.15%), maltose + isomaltose (<0.4%) and maltotriose (<0.2%).

METHOD OF ADMINISTRATION AND APPLICATIONS

The solutions are single-use only and ready to use.

It is advisable to take the preparation cold to facilitate its ingestion.

The method of use depends on the type of test to be performed, whether it is a Glucose Overload test, an O'Sullivan test or a Glucose Stimulation test. Each laboratory will determine how to perform these tests and how the dextrose solution should be used. The dose to be administered may vary depending on the age of the patient. These parameters will be indicated by the medical professional. As an example, for O'Sullivan Test: basal serum glucose should be determined; orally administer a 200mL bottle of Dextrose Solution 50 (orange or lemon); Determine serum glucose at 60 minutes.

CONSERVATION AND STABILITY

Dextrose solutions GLYCOSULL are stable until the expiration date indicated on the label, when kept at room temperature.

Indications of reagent alteration:

These solutions are transparent and free of suspended particles. If the bottle presents turbidity or particles must be discarded.

INTERFERENCES

Glucose is present in most foods and can cause distortions in test result. For this reason they are usually preceded by a period of fasting. Different drugs interferences could alter glucose metabolism affecting the results of the test, so the medical professional should take this into account when evaluating the outcome of such tests, such as thiazide diuretics, oral contraceptives, corticosteroids, beta-blockers, tricyclic antidepressants, steroids, estrogens, thyroxine, heparin, and so on.

PRECAUTIONS

Glucose monohydrate is a widely known substance and it has been widely used in clinical practice for more than 30 years. No toxicological effects have been observed, at suggested doses, with appropriate indications and taking into account the possible contraindications. Contraindications are given by hypersensitivity to any of the components of the product, carbohydrate intolerance, electrolyte imbalance or dehydration states, among others.

Taking these solutions may cause nausea, vomiting, fainting or diarrhea in some patients. For this reason the product must be taken under the supervision of a medical professional.

It is advisable to consult the safety data sheet before handling the product.

The disposal of the residues has to be made according to local regulations.

REFERENCES

Definition, Diagnosis and Classification of Diabetes Mellitus and its Complications. World Health Organization (1999).
Tiezt, NW., Textbook of Clinical Chemistry 5th Edition, W.B. Saunders, Philadelphia (2012). Young, D.S.; Effects of drugs on Clinical Laboratory Tests. 4th Ed. AAC Press. 2000.

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