



Contents

VMG Vitrex Medical Group 4

VITREX lab 7

Vitrex Micro Haematocrit Tubes	8
Vitrex End-to-End Pipettes	10
Vitrex Precalibrated Pipettes	12
Vitrex Capillary & Melting Point Tubes	13
Vitrex Blood Gas Tubes	14
Vitrex Pipette Fixed	16
Vitrex Capillary Dispenser	17
Vitrex Octaheel Baby	18
Vitrex Ultrasafe	20
Vitrex Flex 3	22
Vitrex Lite3	24
Vitrex Steel	26
Vitrex Safety Containers	27
Vitrex Blood Collection System	28
Chiravac™ Blood Collection needles	30
Vitrex Serum Clot Activator Tubes	31
Vitrex Serum Gel and Clot Activator Tubes	32
Vitrex Coagulation Tubes	33
Vitrex EDTA Tubes	34
Vitrex Heparin Tubes	36
Vitrex No Additive Tubes	38
Vitrex Glucose Tubes	39
Vitrex Sedimentation Systems	40
Vitrex Alcohol Swabs	41
Vitrex Tourniquets	42

VITREX diabetes 43

Vitrex Soft Blood Glucose Lancets	44
Vitrex Compact Eject	45
Medoject Fine® Pen Needles	46
Vitrex Safety Pen Needles With Dual Automatic Protective Shields	47
Chirana® Insulin and Tuberculin Syringes	48



CHIRANA injecta

Chirana® 2-part disposable syringes	50
Chirana® 2-part disposable syringes with needles	50
Chirana® 3-part disposable syringes	51
Chirana® 3-part disposable syringes with needles	51
Chirana® Perfusion 3-part syringes	52
Chirana® 3-part disposable Catheter syringes	53
Medoject® Hypodermic needles	54
Medoject® Safety Hypodermic needles	55
Medoject® Blunt fill needles	55
Medoject® Special Hypodermic needles	55
Inox® Ophthalmic needles	55
Chiraflex® IV Cannulas	56
Chiraflex® Safety IV Cannulas	56
Chiraflex® Mandrins	57
Chiraflex® Scalp vein sets	57
Chiraplus™ Infusion set	58
Chirahem™ Transfusion set	59
Chiraplus™ Spikes	60
Chiraplus™ Injection ports	60
Chiraplus™ Stoppers	60
Chiraline™ Extension lines	60
Chiraway™ 3-way stopcock	60
Chiraskin™ Gloves	61
Examination products	61

CHIRMAX surgical

Surgical Meshes	63
Surgical Sutures - Absorbable threads	64
Surgical Sutures - Non-absorbable threads	66





VMG

**Vitrex
Medical
Group**



VMG Vitrex Medical Group

CTI group has per November 2024 been renamed to “VMG Vitrex Medical Group”. The strategy of the group remains the same - to develop and distribute healthcare products of great quality. The vision of the group has changed and instead of being a group of individual companies, the vision is to create “one team” by migrating the companies’ in terms of strategy, values, departments, processes, know-how and capabilities.

As VMG Vitrex Medical Group, we strive to be the leading European medical device & IVD manufacturing group, developing, producing and supplying high-quality healthcare products and tailor-made med-tech solutions that embody special expertise, reliability, trustworthiness and are made with respect and dedication to enhance the life of healthcare professionals and their patients worldwide.

Manufacturing sites

CHIRANA T. Injecta, Slovakia

CHIRANA T. Injecta is a global, state of the art manufacturer of healthcare products, which was founded in 1935.

In 1947 CHIRANA T. Injecta started manufacturing glass syringes and they have ever since focused on development and production of high quality, disposable syringes, needles and metal tubes for the healthcare sector and med-tech companies.

VITREX Medical, Denmark

Vitrex Medical was founded in 1967, where they started producing resistors for televisions, which later on evolved to producing capillary tubes to transport liquids.

With more than 50 years of experience Vitrex Medical has become the largest manufacturer of capillary tubes in the world today by developing equipment, establishing technologies, and optimizing production processes.

CHIRANA T. Injecta, Czech Republic

CHIRANA T. Injecta was founded in 1930. Since the early 1990s they have made substantial investments in upgrading the production resources, enabling the company to offer products in a unique quality.

Their more than 90 years of experience with manufacturing sutures has resulted in top priority within quality control, product development and sutures manufactured in line with state of the art technology.

VITREXlab

Member of Vitrex Medical Group

New and higher standards in the diagnostic laboratory field require an improved performance of laboratory products. Vitrex Lab products and accessories are designed to fulfill the requirements of lab procedures, patient's comfort and user's safety. In addition to the comprehensive, standardized range of products, we have developed a large variety of specialized products.

Range of products:

- Capillary tubes and Haematocrit tubes
- Micro-pipettes
- Safety lancets
- Vacuum tubes and Blood collection needles
- Alcohol swabs and Safety containers

VITREXdiabetes

Member of Vitrex Medical Group

Diabetes is a huge global challenge. There are 415 million people with diabetes worldwide and by 2040 every tenth person is expected to have diabetes. An effective treatment also involves as pleasant a treatment as possible. We aim to deliver quality healthcare solutions for diabetics to meet the demand of both professionals and individual patients. All people are different, so it requires a carefully selected portfolio of products.

Range of products:

- Pen needles
- Dual safety pen needles
- Diabetic lancets
- Insulin syringes

CHIRANAinjecta

Member of Vitrex Medical Group

Drug delivery products are used for medical purposes with the intention of diagnosing, preventing or treating a disease or sickness. These products contribute to help improving the quality of life for patients. We offer a comprehensive range of high-quality drug delivery products for the healthcare sector, which vary in intended use and indications for use.

Range of products:

- 2-part and 3-part syringes
- Perfusion syringes
- Hypodermic needles
- IV administration products
- Examination products

CHIRMAXsurgical

Member of Vitrex Medical Group

For different specialties, procedures and tissue types within human, veterinary and dental market, Chirmax Surgical provides state of the art surgical sutures and meshes to improve the lives of patients worldwide. Our surgical sutures are manufactured in Europe and each individual production process is subject to strict scientific quality control.

Range of products:

- Surgical sutures, ligatures and meshes
- Veterinary bottles and eye needles
- Microsurgery and cardiovascular sutures
- Antibacterial sutures

CHIRANAmed-tech

Member of Vitrex Medical Group

Our CDMO solutions enables companies to leverage our advanced manufacturing capabilities to produce high-quality medical devices under their own brand. We offer versatile capabilities and technologies within glass, plastic, metal and wire processing as well as plastic component molding, assembly, packing, sterilization and finishing.

Examples of CDMO projects:

- Specialty drug delivery cannulas
- Pre-fill syringe cannulas
- Plunger rod for PFS
- Capillary tubes for IVF procedures
- Medical procedure kits for e.g. vaccination





 **VITREX lab**
Member of Vitrex Medical Group



Vitrex Micro Haematocrit Tubes

Vitrex Micro Haematocrit Tubes are quality products, that conform to international standard ISO 12772. The tubes are firepolished at the colour-coded end to reduce wear on the centrifuge seals. We recommend the colour-coded end to be sealed with Vitrex Wax plate. The colour-coding on the tubes and packaging indicates the impregnation of the tube: plain (blue colour), sodium heparin (red colour), ammonium heparin (green colour) and lithium heparin (orange colour). Each heparinised tube contains at least 4-7 IU heparin as prescribed by the standard. Vitrex Micro Haematocrit tubes are also available in plastic with plain, sodium or ammonium heparin in a volume of ~50 µl. The plastic tube provides the same effect as the glass tube.

The tubes are available in two types of packaging: Dispenser vials of 100 tubes each for hygienic and easy handling; and economy cardboard boxes of 1000 tubes, which reduce transport and handling costs. Special requirements regarding marking, dimensions or impregnation are possible.

Accessories

Vitrex Wax plate has 2x24 numbered positions per pair, which makes it easy to seal and transport samples. Each pair can be forced apart and used separately.

Product	Article no.
Wax plate 2 x 5 pcs.	140013





Standard Micro Haematocrit Tubes

Colour Code	Impregnation	d (mm) +/- 0.05 mm	D (mm) +/- 0.05 mm	L (mm) +/- 0.5 mm	Vol. (µl)	Article No. 10 x 100	Article No. 1 x 1000
●	Plain	1.15	1.55	75	~75	160213	163213
●	Sodium Hep.	1.15	1.55	75	~75	161313	163613
●	Ammon. Hep.	1.15	1.55	75	~75	160813	164013
●	Lithium Hep.	1.15	1.55	75	~75	161813	*
●	EDTA	1.15	1.55	75	~75	167313	*

Non-standard Micro Haematocrit Tubes

Shelf life 4 years

Colour Code	Impregnation	d (mm) +/- 0.05 mm	D (mm) +/- 0.05 mm	L (mm) +/- 0.5 mm	Vol. (µl)	Article No. 10 x 100	Article No. 1 x 1000
●	Plain	0.95	1.35	75	~50	160413	163413
●	Sodium Hep.	0.95	1.35	75	~50	161513	163813
●	Ammon. Hep.	0.95	1.35	75	~50	160913	*
●	Lithium Hep.	0.95	1.35	75	~50	160613	*

Precalibrated Micro Haematocrit Tubes

Shelf life 4 years

Colour Code	Impregnation	d (mm)	Tol. d. mm	D (mm) +/- 0.05 mm	L (mm) +/- 0.5 mm	Vol. 1 (µl)	Vol. 2 (µl)	Article No. 10 x 100
●	Sodium Hep.	0.55	+/-0.05	1.55	75	~14	~18	162713
●	Ammon. Hep.	0.55	+/-0.05	1.55	75	~14	~18	162813
●	Lithium Hep.	1.023	+/-0.008	1.50	75	~32	~61	161113


Plastic Micro Haematocrit Tubes

Shelf life 4 years

Colour Code	Impregnation	d (mm) +/- 0.08 mm	D (mm) +/- 0.08 mm	L (mm) ± 1.0mm	Vol. (µl)	Article No. 10 x 100
●	Plain	0.85	1.55	75	~50	260213
●	Sodium Hep.	0.85	1.55	75	~50	261313
●	Ammon. Hep.	0.85	1.55	75	~50	260813

* Available upon request
For special tubes contact us

Shelf life 2 years





Vitrex End-to-End Pipettes

Vitrex End-to-End Pipettes are intended for applications where accurate capillary tube volume is required. The pipettes are in conformity with ISO 7550.

The pipettes are packed in glass vials with a dispenser top, that makes it easy to shake out one tube at a time and thereby ensuring hygienic handling.

The dispenser tops are colour-coded to indicate the tube volume in conformity with ISO 1769. Ten vials are packaged in a sturdy carton in order to protect the tubes from damage during handling or transport.

Vitrex tolerances

Vol. (µl)	Accuracy	CV
1 - 10	± 0.50%	0.75%
20 - 100	± 0.50%	0.60%

Accessories

Product	Article No.	
Pipette holder - Standard	140413	
Pipette holder - Compact	140313	





Technical Specifications

Colour Code	Vol. (µl)	L (mm) +/- 1mm	Pcs./box	Plain	Natrium Hep.	Ammon. Hep.	K2EDTA
●	1	29	10x100	165013	167513 *	169013 *	*
●	2	29	10x100	165113	167613 *	169113 *	*
●	3	29	10x100	165213	167713 *	169213 *	*
●	4	29	10x100	165313	167813 *	169313 *	*
○	5	29	10x100	165413	167913 *	169413 *	170113 *
●	10	29	10x100	174313	177113 *	179213 *	174113 *
●	20	29	10x100	174513	177213 *	179313 *	174413 *
○	25	29	10x100	174613	177313 *	179413 *	178213 *
●	40	29	10x100	174813	177413 *	179513 *	*
●	44.7	29	10x50	175013	177713 *	178613 *	*
●	50	29	10x50	174913	177513 *	179613 *	171113 *

* Available upon request
For lithium hep. please contact us

Shelf life 4 years





Vitrex Precalibrated Pipettes

Vitrex Precalibrated Pipettes are intended for applications where accurate capillary tube volume is required. The pipettes are in conformity with ISO 7550.

The pipettes are packed in cardboard vials. The pipettes have a standard length of 125 mm and the colour-coding on each pipette indicates the volume, which is measured from the base and up to the thin black ringmark.

Vitrex tolerances

Vol. (µl)	Accuracy	CV
1-5	± 0.30%	0.60%
10-200	± 0.25%	0.50%

Technical Specifications

														
Vol. (µl)	1,2,3,4,5	5	10	5-10	20	25	40	20-40	44.7	50	20-50	100	50-100	100-200
Pcs.	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x250	4x125
Plain	128813	126013	126213	125513	126413	126613	129313	125913	126813	127013	125613	127213	125713	129413
Na Hep.	128913	127413	127613	125213	127813	128013	128413	129913	128213	128513	129113	128613	127113	127713

250 pcs. except 100-200 µl, which are packed in vials of 125 pcs.
Other coatings available upon request i.e. Lithium, Ammonium, EDTA.

Shelf life 4 years





Vitrex Capillary and Melting Point Tubes

Vitrex Capillary Tubes are perfect for the handling and transport of small volumes of liquid, melting point determinations and many other purposes. Melting Point tubes are uniform and cut at right angles. The tubes are packed as bulk in boxes. The tubes are also possible to pack in vials of 100 pcs. upon special request.



Technical Specifications

d (mm) +/- 0.05	D (mm) +/- 0.05	L (mm) +/- 1	Pcs.	Capillary Tubs Both ends open	Melting Point One end Closed
0.80	1.00	80	2000	210013	218013
0.80	1.00	100	2000	210113	218113
0.90	1.20	70	1000	210213	218213
0.90	1.20	100	1000	210313	218313
1.00	1.40	75	1000	210413	218413
1.00	1.40	80	1000	210513	218513
1.00	1.40	100	1000	210613	218613
1.00	1.40	125	1000	210713	218713
1.00	1.40	150	1000	210813	218813
1.00	1.40	200	5000	210913	218913
1.15	1.55	80	1000	211013	219013
1.15	1.55	100	1000	211113	219113
1.20	1.65	75	500	211213	219213
1.20	1.65	100	500	211313	219313
1.20	1.65	125	1000	211413	219413
1.20	1.65	150	1000	211513	219513
1.20	1.65	200	3000	211613	219613
1.40	1.90	125	500	211713	219713
1.50	2.00	75	500	211813	219813
1.50	2.00	100	500	211913	219913

For other dimensions, tolerances, sealing, impregnation or packaging, please contact us.

Shelf life 4 years




Vitrex Blood Gas Tubes

Vitrex capillary tubes for blood gas analyses are made from heavy walled glass, suitable for 2.3 mm Ø adaptors. The uniform heparinisation secure proper function of the tubes for blood gas analysis. Vitrex Blood Gas Tubes are also available in plastic with sodium heparin in the same sizes as Vitrex Blood Gas Tubes in glass.

Accessories

Stirrers packed in vials of 250 pcs.
End Caps packed in boxes of 500 pcs.
Magnets packed in bags of 10 pcs.





Blood Gas Tubes

Vol. (µl)	Pcs.	d (mm) +/- 0.05 mm	D (mm) +/- 0.05 mm	L (mm) +/- 1 mm	Na Hep. Article No.	Li Hep. Article No.
130 (min)	4 x 250	1.60	2.30	75	182513	182613
175 (min)	4 x 250	1.60	2.30	100	182413	183313
220 (min)	4 x 250	1.60	2.30	125	182313	182013

Special types upon request.

Shelf life 4 years

Plastic Blood Gas Tubes

Vol. (µl)	Pcs.	d (mm) +/- 0.05 mm	D (mm) +/- 0.05 mm	L (mm) +/- 1 mm	Plastic Na Hep. Article No.
130 (min)	4 x 200	1.60	2.30	75	282513
175 (min)	4 x 200	1.60	2.30	100	282413
220 (min)	4 x 200	1.60	2.30	125	282313

Other types available upon request.

Shelf life 2 years

Accessories

Description	For tubes with D	Pcs.	Article No.
End caps for D (mm)	2.20 - 2.80 mm	10 x 500	140813
End caps for D (mm)	1.50 - 2.20 mm	10 x 500	140713
Stirrers		10 x 250	140913
Magnets		1 x 10	142400





Vitrex Pipette Fixed

LOW COST

This miniature fixed volume pipette is specially designed for use with Diagnostic Kits.

HIGH ACCURACY AND PRECISION

Each pipette is individually calibrated as per the ISO 8655 standard.

COLOUR-CODING

The colour-coding for each pipette indicates the volume - ensuring easy identification in the laboratory environment.

TWO STEP PLUNGER OPERATION

This unique feature makes the pipette a genuine miniature micropipette.

PRACTICAL SIZE

At a length of 130 mm this pipette is the optimum practical size for the pipetting function.

FULLY AUTOCLAVABLE

Autoclavable at 121°C & 15 psi.

Recommended use

Stage 1:

Depress the plunger to the first stage to aspirate the fluid. Place the tip of the micropipette in the desired fluid. Slowly release the pressure and allow the plunger to rise allowing the fluid to aspirate.

Stage 2:

Transfer the fluid to the desired location by again depressing the plunger to the first stage. Press the plunger in full (second stage) in order to fully dispense the fluid.

This operation also allows the Vitrex Pipette Fixed to be used for both reverse and repetitive pipetting.

Technical Specifications

Colour Code	Vol. µl	Accuracy ±%	Precision CV%	Pcs.	Article No.
White	5	2.5	1.5	1	560013
Orange	10	1.2	0.8	1	561013
Black	20	1	0.5	1	562013
Grey	25	1	0.4	1	563013
Green	50	1	0.4	1	564013
Blue	100	0.8	0.3	1	565013
Red	200	0.8	0.3	1	566013




Step 1



Step 2



Step 3



Step 4

Vitrex Capillary Dispenser

Vitrex Capillary Dispenser is a unique liquid handling tool, developed to empty capillary tubes for blood or other liquids.

FULLY AUTOCLAVABLE

Autoclavable at 121°C & 15 psi for sensitive applications.

LOW COST

This dispenser is specially designed for use with Diagnostic Kits.

UNIQUE TIP CONE

Designed to accept most Vitrex Capillary Tubes.

OPTIMUM SIZE

At a length of 130 mm this dispenser is the optimum practical size for blood collection.

Technical Specifications

Colour Code	Vol. µl	Pcs.	Article No.
Grey	Most Vitrex Capillary Tubes	1	569013

Recommended use

Step 1: Loading the device

Insert an unused capillary tube into the pipette tip. The end of the capillary tube must be visible in the tip.

Step 2: Filling the capillary tube

Hold the pipette horizontally and let the capillary tube be filled by capillary effect. Be sure not to press the plunger of the pipette.

Step 3: Emptying the capillary tube

Empty the capillary tube by pressing the plunger gently to the first stage. The liquid will be dispensed, but a small amount of liquid will be left in the tube. After a pause continue to depress the plunger to the second stage. The small amount of liquid left in the tube has also dispensed.

Step 4: After use

Remove the capillary tube carefully and discard in a suitable safety container.



Vitrex Octaheel Baby

Sterile Cutting Safety lancets

Vitrex Octaheel Baby lancets are intended for high volume blood flow with minimum depth of penetration. The safety lancets are primarily applicable for newborns and preemies. However, they can also be recommended for other applications, where short penetration depth and/or minimal compression of tissue is important.

Vitrex Octaheel Baby lancet is the latest edition offering improved ergonomics.



ACTIVATION BUTTON: Secured by safety tab.

NO CHOCKING HAZARD: Non-detachable safety tab.

USER SAFETY: Blade is fully shielded before and after use.

EASIER HANDLING: Anti-slip texture and one hand operation.

USER CONVENIENCE: Pre-loaded, button activated design.

PATIENT SAFETY: Single-use only. Sterilised by gamma irradiation.

PATIENT COMFORT: Ultra sharp blade with high speed cutting securing minimal trauma and fast healing.

NARROWED LANCING SLOT: Easier positioning.



Recommended use

1. Release the trigger lock mechanism by pushing down the safety tab.
2. Reference position for infant heel blood sampling.
3. Raise the baby's foot to a suitable position, place the blade slot flush against the heel site and activate the lancet.
4. Collect the blood sample with a capillary tube.
5. Dispose the safety lancet in a suitable safety container.

Technical Specifications

Colour Code		Specification	Article No.	Variant	Indicative use
Green		0.65 mm Depth / 1.40 mm Incision Length	335313	Micro-Preemie	Typically refers to premature infants with a birth weight below 1000 grams
Orange		0.85 mm Depth / 1.75 mm Incision Length	334313	Preemie	Newborns with a birth weight below 2500 grams, most of whom are premature
Yellow		1.00 mm Depth / 2.50 mm Incision Length	333313	Newborn	Newborns with a birth weight between 2500 grams and 4000 grams
Pink		1.14 mm Depth / 2.80 mm Incision Length	332313	Newborn Plus	Generally refers to newborns with a birth weight over 4000 grams
Blue		2.00 mm Depth / 3.00 mm Incision Length	331313	Toddler	Refers to children aged 1 to 3 years old, and does not pertain to birth weight



Transit packaging

Pallet

Packaging:

Blister Packed

Primary packaging: 50 pcs.
 Transit packaging: 200 pcs.
 Pallet: 23,200 pcs.

Shelf life 5 years



Vitrex Ultrasafe

Sterile Safety lancet

Vitrex Ultrasafe lancet is intended for blood glucose and other in vitro diagnostic tests. Health workers worldwide face accidental needle-stick injury while caring for patients, putting themselves at risk to bloodborne illnesses including the AIDS virus and Hepatitis.

PRESSURE ACTIVATED

UNIQUE DESIGN: Ergonomically optimized for better handling

HIGH QUALITY SILICONE COATED NEEDLE: Less pain and smooth penetration

EASY TARGETING OF PUNCTURE SITE




Unique design, ergonomically optimized for better handling



1



2



3

Recommended use

1. Carefully twist off the protective cap.
2. Select puncture site and press the device firmly against the site to release the lancet.
3. Dispose the safety lancet in a suitable safety container.
4. Collect the blood sample with a capillary tube.

Technical Specifications

Colour Code	Specification	Penetration Depth	Article No.	Indicative use
Green	17G blade	2.0 mm	393413	For heel puncture in average-weight newborns and older children
Blue	21G needle	1.8 mm	397413	Diabetes clinics/hospital wards for blood glucose monitoring
Orange	21G needle	2.2 mm	395413	Haemoglobin (blood banks), cholesterol, blood typing, coagulation, blood gas and other diagnostic tests
Pink	21G needle	2.8 mm	396413	As the Orange lancet, but for thicker adult skin
Yellow	26G needle	1.8 mm	394413	As the Orange lancet, but lower blood flow and less pain
Rose	28G needle	1.5 mm	392413	Diabetes clinics/hospital wards for blood glucose monitoring
Purple	28G needle	1.8 mm	398413	Low blood flow, less pain
Turquoise	30G needle	1.8 mm	399413	Diabetes clinics/hospital wards for blood glucose monitoring. For minimal pain

Shelf life 5 years



Transit packaging

Pallet

Packaging:

Primary packaging:

100 pcs.

Transit packaging:

1,000 pcs.

Pallet:

98,000 pcs.



Vitrex Flex 3

Depth Adjustable Safety lancet

Vitrex Flex 3 safety lancet is an All-in-One lancet for emergency use due to its adjustable design, which includes three different penetration depths for a single unit. The lancet is intended for blood glucose and other in vitro diagnostic tests.

USER CONVENIENCE: Pre-loaded and pressure activated design. The trigger position indicates the before and after use.

USER SAFETY: Needle is fully shielded before and after use.

COMFORT: High quality tri-bevel needle with high speed penetration.

PATIENT SAFETY: Single-use only. Sterilised by gamma irradiation.



Avoid misfiring when adjusting the depth. Ergonomic design for improved holding.



Improved safety cap removal
– cap automatically semi
ejects, when turned
to remove.

Technical Specifications

Colour Code	Specification	Penetration Depth	Article No.
Blue	23G	1.30 mm - 1.80 mm - 2.30 mm	395813
Yellow	26G	1.30 mm - 1.80 mm - 2.30 mm	396813
Orange	28G	1.30 mm - 1.80 mm - 2.30 mm	397813

Shelf life 5 years




Transit packaging

Pallet

Packaging:

Primary packaging:
Transit packaging:
Pallet:

100 pcs.
1,000 pcs.
84,000 pcs.



Penetration depth adjustment
is done without touching the
activation mechanism.



Vitrex Lite3

Sterile Safety lancet

Vitrex Lite3 lancet is intended for blood glucose and other in vitro diagnostic tests. Health workers worldwide face accidental needle-stick injury while caring for patients, putting themselves at risk to bloodborne illnesses including the AIDS virus and Hepatitis.

BUTTON ACTIVATION

USER SAFETY: Needle is fully shielded before and after use.


USER CONVENIENCE: Pre-loaded, pressure activated design.

PATIENT COMFORT: High quality, ultra sharp tri-bevel needle or bi-bevel blade, high speed penetration.

EFFICACY: Variety of point configurations and penetration depths to suit most capillary blood requirements.

PATIENT SAFETY: Single-use only. Sterilised by gamma irradiation.





Recommended use

1. Twist off the protective cap.
2. Select fingertip insertion site slightly off centre to minimise any discomfort and maximise blood flow. Hold the lancet firmly against site and press the release button to activate the device.
3. Gently massage the finger in the direction of the puncture site to obtain the required volume of blood.
4. Dispose the safety lancet in a suitable safety container.

Technical Specifications

Colour Code	Specification	Penetration Depth	Article No.	Indicative use
Green	18G blade	1.8 mm	373213	For heel puncture in average-weight newborns and older children
Orange	21G needle	1.8 mm	377213	Haemoglobin (blood banks), cholesterol, blood typing, coagulation, blood gas and other diagnostic tests
Pink	21G needle	2.4 mm	378213	As the Orange lancet, but for thicker skin
Dark Blue	23G needle	1.8 mm	379213	Blood glucose testing and tests requiring larger samples
Yellow	26G needle	1.8 mm	375213	Diabetes clinics/hospital wards for blood glucose monitoring
Blue	26G needle	2.4 mm	376213	As the Yellow lancet, but for thicker skin
Rose	28G needle	1.5 mm	371213	Diabetes clinics/hospital wards for blood glucose monitoring
Purple	28G needle	1.8 mm	374213	As the Rose lancet, but for thicker skin
Beige	30G needle	1.8 mm	372213	Diabetes clinics/hospital wards for blood glucose monitoring. For minimal pain

Shelf life 5 years



Transit packaging

Pallet

Packaging:

Primary packaging: 100 pcs.
Transit packaging: 1,000 pcs.
Pallet: 90,000 pcs.



Vitrex Steel

General Purpose Capillary Blood Lancet

Vitrex Steel lancets are traditional handheld, flat and stainless steel blood lancets for general purpose use. The steel lancets are available in standard and in allergy version.

LOW COST

Simple, flat, steel lancet.

SAMPLE VOLUME

Excellent blood flow.

HIGH QUALITY

Hardened surgical steel.

MINIMUM DISCOMFORT

Modern asymmetric design.

STERILE

Single-use only. Sterilised by gamma irradiation.

Technical Specifications

Product	Description	Article No.
Vitrex Steel	Flat blade sterile lancet	350213
Vitrex Steel - allergy	Flat blade sterile lancet	365213

Shelf life 5 years



Packaging:
 Primary packaging: 200 pcs.
 Secondary packaging: 10,000 pcs.
 Transit packaging: 40,000 pcs.
 Pallet: 720,000 pcs.



Vitrex Safety Containers

Vitrex Safety Containers are designed for hospitals and healthcare institutions for clinical disposal of lancets, needles, scalpels, syringes and other types of sharps.

Each type of lid is specially designed to help the user dispose needles, syringes, etc. lid with reversible closure can be safely fixed.

Technical Specifications

Description	Size	Content	Article No.
Safety container	4.7 cm x 12.9 cm	0.2 l.	193800
Safety container	9 cm x 10.8 cm	0.5 l.	194000
Safety container	13.7 cm x 12.5 cm	1.5 l.	194200
Safety container	13.7 cm x 18 cm	2 l.	194300
Safety container	21.1 cm x 19.5 cm	5 l.	194700



4.
Blood
Collection
Needle

3.
Blood
Collection
Needle
Holder

2.
Safety
Cap

1.
Blood
Collection
Tube



Vitrex Blood Collection System

New and higher standards in the diagnostic laboratory field require an improved performance of the blood collection system. Vitrex Blood Collection System is the answer.

Vitrex Blood Collection System consists of a full range of products including our Blood Collection Tubes, Needles and Holders, as well as the newly added products: Luer Slip adapter, Luer Slip adapter + holder, Luer Slip adapter + Intravenous needle and Luer Slip adapter + holder + Intravenous needle.

Vitrex offers two different brands of blood collection tubes: Vacusence and Vacusence+. Vacusence tubes are of Turkish origin with a lead time of 8-10 weeks. Vacusence+ tubes are made in Europe with a lead time of 3-4 weeks. Vitrex Blood Collection Tubes are made of unbreakable PET tube, butyl rubber stopper and polyethylene cap. We primarily offer medical paper labelled tubes, but tubes with transparent label or prebarcoded tubes are also available upon request.

Letter Code	Additive
K2E	EDTA Dipotassium Salt
K3E	EDTA Tripotassium Salt
9NC	Trisodium Citrate 9:1
4NC	Trisodium Citrate 4:1
FE	Sodium Fluoride I EDTA
FX	Sodium Fluoride/Potassium Oxalate
LH	Lithium Heparin
NH	Sodium Heparin
CAT	Clot Activator
Z	None



All components of Vitrex Blood Collection System are designed to fulfill the requirements of blood collection procedures, improve patient's comfort, help user's safety and get reliable test results. The latest products in our Blood Collection System range are of corresponding high quality and carefully selected to complete the range.

USER VERSATILITY: Series combination blood collection types for customer selection.

USER CONVENIENCE: Allows for multiple draws with a single vein puncture.

EXCELLENT PERFORMANCE: High quality rubber sleeve has excellent self-sealing performance.

FLEXIBLE PACKAGING: Single blister package, single pouch or single PE bag.

Description	Colour Code	Size	Diameter (mm)	Pcs./ Box	Pcs./ Carton	Article No.	Article No. Safety needle
Luer Slip adapter + Intravenous needle	Green	21G	0.80 x 19mm	100	2000	VBCS72221	VBCS73221S
Luer Slip adapter + Intravenous needle	Blue	23G	0.60 x 19mm	100	2000	VBCS72223 *	VBCS73223S *
Luer Slip adapter + holder + Intravenous needle	Green	21G	0.80 x 19mm	50	1000	VBCS82221	VBCS83221S
Luer Slip adapter + holder + Intravenous needle	Blue	23G	0.60 x 19mm	50	1000	VBCS82223 *	VBCS83223S *

* Contact Sales Department for more information.

Shelf life 5 years



Chiravac™ Blood Collection needles

Chiravac Blood Collection needles are made of polypropylene hub and stainless steel tube. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. The needles are intended for blood collection with vacuum tubes. They are individually packed in boxes of 100 pcs, cartons of 1000 pcs. Upon request the needles are also available as safety or invisible version.

USER CONVENIENCE: Colour-coded needle hub, screw-on needle hub compatible with all available holders, easy-to-open packaging and thin wall tube material.

COMFORT: Grinded, precision lancet point and silicone treated tip for smooth penetration. Comfortable collection experience.

PATIENT SAFETY: Single-use only. Sterilised by ETO.

FDA CLEARED

ACCESSORIES

Blood Collection Needle Holders with and without safety function are compatible with our Chiravac Blood Collection Needles and our different Blood Collection Tubes. The holders are packed in PE bags of 50 or 100 pcs.

Colour Code	Size Metric	Size G/"	Article No.
Blue	0.60 x 38 mm	23G x 1 1/2"	CHBCN23112
Black	0.70 x 38 mm	22G x 1 1/2"	CHBCN22112
Green	0.80 x 38 mm	21G x 1 1/2"	CHBCN21112
Yellow	0.90 x 38 mm	20G x 1 1/2"	CHBCN20112
Pink	1.20 x 38 mm	18G x 1 1/2"	CHBCN18112
Blue	0.60 x 25 mm	23G x 1"	CHBCN23100
Black	0.70 x 25 mm	22G x 1"	CHBCN22100
Green	0.80 x 25 mm	21G x 1"	CHBCN21100
Yellow	0.90 x 25 mm	20G x 1"	CHBCN20100
Pink	1.20 x 25 mm	18G x 1"	CHBCN18100

Shelf life 5 years

Description	Pcs./Box	Pcs./Carton	Article No.
Blood Collection Needle Holder	100	1,000	CHBCNH01
Blood Collection Needle Safety Holder	50	800	CHBCNHS01

Shelf life 5 years



Vitrex Serum Clot Activator Tubes

Vitrex Serum Clot Activator Tubes are used for clinical chemistry, immunology, protein electrophoresis, serology, microbiology and toxicology. The inner wall of the Serum Clot Activator tube is specially coated with clot activator, which activates the coagulation process.



Transit packaging

Pallet

Packaging:

Primary packaging: 50/100 pcs.
 Transit packaging: 800/1,000 pcs.
 * Pallet: 54,000 pcs.
 ** Pallet: 48,000 pcs.
 *** Pallet: 38,400 pcs.

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Article No.
Serum CAT/Clot Activator	1	13 x 75	Serum Clot Activator	Paper	101201+ *
Serum CAT/Clot Activator	2	13 x 75	Serum Clot Activator	Paper	102202+ *
Serum CAT/Clot Activator	3	13 x 75	Serum Clot Activator	Paper	103203+ *
Serum CAT/Clot Activator	3	13 x 100	Serum Clot Activator	Paper	108203+ **
Serum CAT/Clot Activator	4	13 x 75	Serum Clot Activator	Paper	104204+ *
Serum CAT/Clot Activator	4	13 x 100	Serum Clot Activator	Paper	108204+ **
Serum CAT/Clot Activator	5	13 x 75	Serum Clot Activator	Paper	109205+ *
Serum CAT/Clot Activator	5	13 x 100	Serum Clot Activator	Paper	108205+ **
Serum CAT/Clot Activator	6	13 x 100	Serum Clot Activator	Paper	105206+ **
Serum CAT/Clot Activator	7	13 x 100	Serum Clot Activator	Paper	108207+ **
Serum CAT/Clot Activator	8	16 x 100	Serum Clot Activator	Paper	108208+ ***
Serum CAT/Clot Activator	9	16 x 100	Serum Clot Activator	Paper	106209+ ***
Serum CAT/Clot Activator	10	16 x 100	Serum Clot Activator	Paper	107200+ ***

Shelf life 18 months



Vitrex Serum Gel and Clot Activator Tubes

Vitrex Vacusence Serum Gel and Clot Activator Tubes contain inert gel; during centrifugation, the gel forms a stable barrier between serum and blood cells. The tubes allow to obtain more serum than regular tubes due to the cleaner separation of serum from the clot. The range of application is: clinical chemistry, immunology, protein electrophoresis, serology, microbiology and toxicology.



Transit packaging	
Packaging:	
Primary packaging:	50/100 pcs.
Transit packaging:	800/1,000 pcs.
* Pallet:	54,000 pcs.
** Pallet:	48,000 pcs.
*** Pallet:	38,400 pcs.

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Article No.
Serum Gel and CAT/Clot Activator	1	13 x 75	Serum Gel and Clot Activator	Paper	102301+ *
Serum Gel and CAT/Clot Activator	2	13 x 75	Serum Gel and Clot Activator	Paper	101302+ *
Serum Gel and CAT/Clot Activator	3	13 x 75	Serum Gel and Clot Activator	Paper	102303+ *
Serum Gel and CAT/Clot Activator	3	13 x 100	Serum Gel and Clot Activator	Paper	102330+ **
Serum Gel and CAT/Clot Activator	3.5	13 x 75	Serum Gel and Clot Activator	Paper	103303+ *
Serum Gel and CAT/Clot Activator	4	13 x 75	Serum Gel and Clot Activator	Paper	106304+ *
Serum Gel and CAT/Clot Activator	4	13 x 100	Serum Gel and Clot Activator	Paper	102304+ **
Serum Gel and CAT/Clot Activator	5	13 x 100	Serum Gel and Clot Activator	Paper	104305+ **
Serum Gel and CAT/Clot Activator	6	13 x 100	Serum Gel and Clot Activator	Paper	102306+ **
Serum Gel and CAT/Clot Activator	8	16 x 100	Serum Gel and Clot Activator	Paper	105308+ ***
Serum Gel and CAT/Clot Activator	8.5	16 x 100	Serum Gel and Clot Activator	Paper	102308+ ***

Shelf life 18 months



Vitrex Coagulation Tubes

Vitrex Coagulation Tubes are used to collect venous blood for coagulation tests. Coagulation Tubes have a buffered sodium citrate solution of 0.109 mol/l (3.2%) or 0.129 mol/l (3.8%). Double Tube technology provides an excellent condition for the test of PT and APTT values. The range of application is: Hemostasis system diagnostics.



Packaging:
 Primary packaging: 50/100 pcs.
 Transit packaging: 800/1,000 pcs.
 * Pallet: 54,000 pcs.
 ** Pallet: 48,000 pcs.
 *** Pallet: 38,400 pcs.

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	3.2%	3.8%
					Article No.	Article No.
Coagulation 9NC/Sodium Citrate	0.9	13 x 75	Sodium Citrate	Paper	101801+ *	106801+ *
Coagulation 9NC/Sodium Citrate	1.8	13 x 75	Sodium Citrate	Paper	102802+ *	106802+ *
Coagulation 9NC/Sodium Citrate, Double Wall	1.8	13 x 75	Sodium Citrate	Paper	1028020+ *	1068020+ *
Coagulation 9NC/Sodium Citrate	1.8	13 x 100	Sodium Citrate	Paper	103802+ **	106820+ **
Coagulation 9NC/Sodium Citrate	2.7	13 x 75	Sodium Citrate	Paper	103803+ *	107803+ *
Coagulation 9NC/Sodium Citrate, Double Wall	2.7	13 x 75	Sodium Citrate	Paper	1038030+ *	1078030+ *
Coagulation 9NC/Sodium Citrate	2.7	13 x 100	Sodium Citrate	Paper	103830+ **	107830+ **
Coagulation 9NC/Sodium Citrate	3.6	13 x 75	Sodium Citrate	Paper	103840+ *	108804+ *
Coagulation 9NC/Sodium Citrate	3.6	13v100	Sodium Citrate	Paper	103804+ **	108840+ **
Coagulation 9NC/Sodium Citrate	4.5	13 x 75	Sodium Citrate	Paper	103850+ *	108805+ *
Coagulation 9NC/Sodium Citrate	4.5	13 x 100	Sodium Citrate	Paper	103805+ **	108850+ **
Coagulation 9NC/Sodium Citrate #	8.1	16 x 100	Sodium Citrate	Paper	103808+ ***	108808+ ***

Shelf life 9 months
 # 8.1 ml: Shelf life 12 months
 Double Wall: Shelf life 13 months



Vitrex EDTA Tubes

Vitrex EDTA Tubes are used for whole blood analysis within haematology. The range of application is: hematology, PCR analysis, blood typing. EDTA and Gel is mainly used for molecular diagnostics. The gel forms a resistant barrier between the cellular components of blood and serum during centrifugation.

K2EDTA					Plain	+ Gel
Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Article No.	Article No.
K2E/K2EDTA	1	13 x 75	K2EDTA	Paper	101901+ *	112201+ *
K2E/K2EDTA	2	13 x 75	K2EDTA	Paper	102902+ *	112202+ *
K2E/K2EDTA	3	13 x 75	K2EDTA	Paper	103903+ *	111203+ *
K2E/K2EDTA	3	13 x 100	K2EDTA	Paper	103930+ **	113203+ **
K2E/K2EDTA	3.5	13 x 75	K2EDTA	Paper		112235+ *
K2E/K2EDTA	4	13 x 75	K2EDTA	Paper	104904+ *	112240+ *
K2E/K2EDTA	4	13 x 100	K2EDTA	Paper	105904+ **	112204+ **
K2E/K2EDTA	5	13 x 75	K2EDTA	Paper	105905+ *	
K2E/K2EDTA	5	13 x 100	K2EDTA	Paper	105950+ **	112205+ **
K2E/K2EDTA	6	13 x 100	K2EDTA	Paper	106906+ **	112206+ **
K2E/K2EDTA	7	13 x 100	K2EDTA	Paper	106907+ **	
K2E/K2EDTA	8	16 x 100	K2EDTA	Paper	106908+ ***	113208+ ***
K2E/K2EDTA	8.5	16 x 100	K2EDTA	Paper		113209+ ***
K2E/K2EDTA	9	16 x 100	K2EDTA	Paper	107909+ ***	113290+ ***
K2E/K2EDTA	9.5	16 x 100	K2EDTA	Paper	108909+ ***	

Shelf life 18 months



Transit packaging

Pallet

Packaging:
 Primary packaging: 50/100 pcs.
 Transit packaging: 800/1,000 pcs.
 * Pallet: 54,000 pcs.
 ** Pallet: 48,000 pcs.
 *** Pallet: 38,400 pcs.

K3EDTA

Plain + Gel

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Article No.	Article No.
K3E/K3EDTA	1	13 x 75	K3EDTA	Paper	111101+ *	122101+ *
K3E/K3EDTA	2	13 x 75	K3EDTA	Paper	112102+ *	122102+ *
K3E/K3EDTA	3	13 x 75	K3EDTA	Paper	113103+ *	123103+ *
K3E/K3EDTA	3	13 x 100	K3EDTA	Paper	113130+ **	123130+ **
K3E/K3EDTA	3.5	13 x 75	K3EDTA	Paper		124135+ *
K3E/K3EDTA	4	13 x 75	K3EDTA	Paper	113104+ *	125140+ *
K3E/K3EDTA	4	13 x 100	K3EDTA	Paper	114104+ **	125104+ **
K3E/K3EDTA	5	13 x 75	K3EDTA	Paper	114105+ *	
K3E/K3EDTA	5	13 x 100	K3EDTA	Paper	114150+ **	125150+ **
K3E/K3EDTA	6	13 x 100	K3EDTA	Paper	116106+ **	126106+ **
K3E/K3EDTA	7	13 x 100	K3EDTA	Paper	116107+ **	
K3E/K3EDTA	8	16 x 100	K3EDTA	Paper	114108+ ***	127108+ ***
K3E/K3EDTA	8.5	16 x 100	K3EDTA	Paper	114185+ ***	128185+ ***
K3E/K3EDTA	9	16 x 100	K3EDTA	Paper	117109+ ***	129109+ ***

Shelf life 18 months



Vitrex Heparin Tubes

Vitrex Heparin Tubes are coated to inhibit clotting. Heparin Tubes are used to produce plasma for biochemical testing. The Heparin Tubes are also available with Gel separator. The range of application is: clinical chemistry and trace mineral testing.

Lithium Heparin

Plain + Gel

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Article No.	Article No.
LH/Lithium Heparin	1	13 x 75	Lithium Heparin	Paper	105401+ *	102701+ *
LH/Lithium Heparin	2	13 x 75	Lithium Heparin	Paper	101402+ *	102702+ *
LH/Lithium Heparin	3	13 x 75	Lithium Heparin	Paper	101403+ *	101703+ *
LH/Lithium Heparin	3	13 x 100	Lithium Heparin	Paper	101404+ **	101730+ **
LH/Lithium Heparin	4	13 x 75	Lithium Heparin	Paper	105404+ *	
LH/Lithium Heparin	4	13 x 100	Lithium Heparin	Paper	102440+ **	102704+ **
LH/Lithium Heparin	5	13 x 75	Lithium Heparin	Paper	105450+ *	
LH/Lithium Heparin	5	13 x 100	Lithium Heparin	Paper	105405+ **	102705+ **
LH/Lithium Heparin	6	13 x 100	Lithium Heparin	Paper	103406+ **	
LH/Lithium Heparin	7	13 x 100	Lithium Heparin	Paper	103407+ **	
LH/Lithium Heparin	8	16 x 100	Lithium Heparin	Paper	105408+ ***	103706+ ***
LH/Lithium Heparin	9	16 x 100	Lithium Heparin	Paper	104409+ ***	102709+ ***

Shelf life 18 months



Transit packaging

Pallet

Packaging:

Primary packaging: 50/100 pcs.

Transit packaging: 800/1,000 pcs.

* Pallet: 54,000 pcs.

** Pallet: 48,000 pcs.

*** Pallet: 38,400 pcs.

Sodium Heparin

Plain

+ Gel

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Article No.	Article No.
NH/Sodium Heparin	1	13 x 75	Sodium Heparin	Paper	104501+ *	104601+ *
NH/Sodium Heparin	2	13 x 75	Sodium Heparin	Paper	104502+ *	104602+ *
NH/Sodium Heparin	3	13 x 75	Sodium Heparin	Paper	104530+ *	101603+ *
NH/Sodium Heparin	3	13 x 100	Sodium Heparin	Paper	104503+ **	104630+ **
NH/Sodium Heparin	4	13 x 75	Sodium Heparin	Paper	101504+ *	
NH/Sodium Heparin	4	13 x 100	Sodium Heparin	Paper	104504+ **	104604+ **
NH/Sodium Heparin	5	13 x 75	Sodium Heparin	Paper	104550+ *	
NH/Sodium Heparin	5	13 x 100	Sodium Heparin	Paper	104505+ **	102605+ **
NH/Sodium Heparin	6	13 x 100	Sodium Heparin	Paper	102506+ **	
NH/Sodium Heparin	7	13 x 100	Sodium Heparin	Paper	102507+ **	
NH/Sodium Heparin	8	16 x 100	Sodium Heparin	Paper	104509+ ***	103608+ ***
NH/Sodium Heparin	9	16 x 100	Sodium Heparin	Paper	103509+ ***	104609+ ***

Shelf life 18 months



Vitrex No Additive Tubes

Vitrex No Additive Tubes do not contain any additives. The tubes are used to collect samples of body fluids (blood, urine, cerebrospinal fluid, exudates, etc.). The tubes can be used for specimen analysis, storage and transport. The range of application is: analysis of body fluids in clinical biochemistry, serology and immunology.



Transit packaging

Pallet

Packaging:
 Primary packaging: 50/100 pcs.
 Transit packaging: 800/1,000 pcs.
 * Pallet: 54,000 pcs.
 ** Pallet: 48,000 pcs.
 *** Pallet: 38,400 pcs.

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Article No.
Z/No Additive	1	13 x 75	No additive	Paper	105101+ *
Z/No Additive	2	13 x 75	No additive	Paper	101102+ *
Z/No Additive	3	13 x 75	No additive	Paper	105103+ *
Z/No Additive	3	13 x 100	No additive	Paper	105130+ **
Z/No Additive	4	13 x 75	No additive	Paper	102204+ *
Z/No Additive	4	13 x 100	No additive	Paper	105104+ **
Z/No Additive	5	13 x 75	No additive	Paper	105150+ *
Z/No Additive	5	13 x 100	No additive	Paper	105105+ **
Z/No Additive	6	13 x 100	No additive	Paper	103106+ **
Z/No Additive	7	13 x 100	No additive	Paper	105107+ **
Z/No Additive	8	16 x 100	No additive	Paper	105108+ ***
Z/No Additive	9	16 x 100	No additive	Paper	104109+ ***
Z/No Additive	10	16 x 100	No additive	Paper	105110+ ***

Shelf life 18 months



Vitrex Glucose Tubes

Vitrex Glucose Tubes contain Sodium Fluoride and Potassium Oxalate. Glucose tubes are used for glucose determinations like blood sugar and sugar tolerance. The tubes contain a glucose stabiliser, which allows you to block the process of glycolysis for up to 24 hours at room temperature. The range of application is: research of glucose and lactate.



Transit packaging Pallet

Packaging:
 Primary packaging: 50/100 pcs.
 Transit packaging: 800/1,000 pcs.
 * Pallet: 54,000 pcs.
 ** Pallet: 48,000 pcs.
 *** Pallet: 38,400 pcs.

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Article No.
FX/ Sodium Flouride + Potassium Oxalate	1	13 x 75	Sodium Flouride + Potassium Oxalate	Paper	111401+ *
FX/ Sodium Flouride + Potassium Oxalate	2	13 x 75	Sodium Flouride + Potassium Oxalate	Paper	111402+ *
FX/ Sodium Flouride + Potassium Oxalate	2	13 x 100	Sodium Flouride + Potassium Oxalate	Paper	113420+ **
FX/ Sodium Flouride + Potassium Oxalate	3	13 x 75	Sodium Flouride + Potassium Oxalate	Paper	113430+ *
FX/ Sodium Flouride + Potassium Oxalate	3	13 x 100	Sodium Flouride + Potassium Oxalate	Paper	113403+ **
FX/ Sodium Flouride + Potassium Oxalate	4	13 x 75	Sodium Flouride + Potassium Oxalate	Paper	112404+ *
FX/ Sodium Flouride + Potassium Oxalate	4	13 x 100	Sodium Flouride + Potassium Oxalate	Paper	113404+ **
FX/ Sodium Flouride + Potassium Oxalate	5	13 x 75	Sodium Flouride + Potassium Oxalate	Paper	113450+ *
FX/ Sodium Flouride + Potassium Oxalate	5	13 x 100	Sodium Flouride + Potassium Oxalate	Paper	113405+ **
FX/ Sodium Flouride + Potassium Oxalate	6	13 x 100	Sodium Flouride + Potassium Oxalate	Paper	113406+ **
FX/ Sodium Flouride + Potassium Oxalate	7	13 x 100	Sodium Flouride + Potassium Oxalate	Paper	113407+ **
FX/ Sodium Flouride + Potassium Oxalate	8	16 x 100	Sodium Flouride + Potassium Oxalate	Paper	113408+ ***
FX/ Sodium Flouride + Potassium Oxalate	9	16 x 100	Sodium Flouride + Potassium Oxalate	Paper	113409+ ***

Shelf life 18 months



Vitrex Sedimentation Systems

Vitrex Sedimentation Tubes contain 3.2% sodium citrate for ESR test. The tubes are available with an open or closed sedimentation system.



Transit packaging

Packaging:
 Primary packaging: 50/100 pcs.
 Transit packaging: 800/1,000 pcs.
 * Pallet: 54,000 pcs.
 ** Pallet: 48,000 pcs.

Pallet

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Article No.
ESR Tube	1.6	13 x 75	4NC/Sodium Citrate 3.2%	Paper	112602+ *
ESR Tube	2	13 x 75	4NC/Sodium Citrate 3.2%	Paper	111602+ *
ESR Tube	2.4	13 x 75	4NC/Sodium Citrate 3.2%	Paper	113624+ *
ESR Tube	2.4	13 x 100	4NC/Sodium Citrate 3.2%	Paper	114624+ **
ESR Tube	3	13 x 75	4NC/Sodium Citrate 3.2%	Paper	114603+ *

Shelf life 9 months



Accessories

Vitrex Accessories are necessary products for blood collection.

Vitrex Alcohol Swabs

Vitrex Alcohol Swabs are an easy and convenient way to cleanse the skin and help fight bacteria at the site of application. Use as an antiseptic skin cleanser before injections and blood sampling.

The size of the Disinfection Swab is 30x60 mm. Each Alcohol Swab is individually sealed and disposable for convenience and safety.

Technical Specifications

Description	Size	Content	Identification No.	Pcs./Box	Article No.
Alcohol Swabs	30 x 60 mm	82% Ethanol	1U8K-YGJD-J103-YQRW	100	521413
Alcohol Swabs	30 x 60 mm	82% Ethanol + 0.5% Chlorhexidine	NQ1P-3GSC-U100-G3GE	100	522413

Shelf life 5 years

Alcohol Swabs are currently approved for sales in: Sweden, Finland, UK, Germany, France and Austria.

(Updated: 2023-10. Please consult our website: vitrexmg.com for new updates in the list of approved countries).

Regulated according to: REGULATION (EU) No 528/2012 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 22 May 2012 concerning the making available on the market and use of biocidal products.



Transit packaging

Pallet

Packaging:
 Primary packaging: 100 pcs.
 Transit packaging: 10,000 pcs.
 Pallet: 400,000 pcs.



Accessories

Vitrex Accessories are necessary products for blood collection.

Vitrex Tourniquets

Vitrex Tourniquets are used for tying the vein during blood collection. They are safe and easy to use. The tourniquets are available in automatic version for both adults and children, but also in a manual, disposable version as one-size. Latex free.

Description	Pcs./Box	Article No.
Vitrex Tourniquets with buckle, Adult	50	TM1150
Vitrex Tourniquets with buckle, Children	50	TM1125
Vitrex Disposable Tourniquets	25 pcs./roll	TM1175

Tourniquets with buckle



Transit packaging

Pallet

Packaging:
 Primary packaging: 50 pcs.
 Transit packaging: 500 pcs.
 Pallet: 9,000 pcs.

Tourniquets Disposable



Transit packaging

Pallet

Packaging:
 Primary packaging: 1 roll
 Transit packaging: 100 rolls
 Pallet: 3,000 rolls



 **VITREX** diabetes
Member of Vitrex Medical Group



Vitrex Soft Blood Glucose Lancets

Vitrex Soft Blood Glucose Lancets are intended for blood glucose measurement. The soft lancets are universal and compatible with most lancing devices especially Vitrex Compact Eject. The range includes standard soft lancets and siliconized soft lancets available in different sizes.

GOOD BLOOD FLOW

Triple ground needle point.

VERSATILE

Fit most commonly used lancing devices including the Vitrex Compact Eject lancing device.

PATIENT SAFETY:

Single-use only. Sterilised by gamma irradiation.

CONVENIENT

Boxes have dispensing feature.

NEW SILICONIZED NEEDLE

For less pain and smoother penetration.

Technical Specifications

Colour Code		Specification	Pcs./Box	Article No.
Light blue		23G (0.64 mm)	200	351213
Red		28G (0.36 mm)	200	352213
Green		30G (0.30 mm)	200	354213
Dark blue		32G (0.26 mm)	200	357213

Shelf life 5 years

Technical Specifications Siliconized Needles

Colour Code		Specification	Pcs./Box	Article No.
Red		28G (0.36 mm)	200	352413
Green		30G (0.30 mm)	200	354413


Shelf life 5 years



Transit packaging

Pallet

Packaging:
 Primary packaging: 200 pcs.
 Secondary packaging: 2,000 pcs.
 Transit packaging: 20,000 pcs.
 Pallet: 480,000 pcs.



Recommended use

1. Variable penetration

Dial the appropriate depth on the device between 1 (min.) and 5 (max.).

2. Simple and quick activation and reactivation

Pull back the grey slider to activate the mechanism.

3. Instant operation

Depress the grey button to operate the device.

4. Lancet removal

Simple and reliable lancet ejection.

Vitrex Compact Eject

Vitrex Compact Eject is a lancing device with a push button ejection. The device is clear adjustable with clear fixed end caps designed for Alternate Site Testing (AST).

VERSATILE

Suitable for use with Vitrex Soft Lancets and other most commonly used lancet brands.

VARIABLE PENETRATION

Nine depth settings from 1 to 5.

CONVENIENT

Simple and quick activation, reactivation and non-touch lancet removal.

Technical Specifications

Product	Description	Packaging	Article No.
Vitrex Compact Eject	9-step adjustable depth	Clear plastic sleeve only	369613

Shelf life 5 years





Medoject Fine® Pen Needles

Medoject Fine Pen Needles are made of polypropylene hub and stainless steel tube. They are intended for Insulin application via most Insulin pens. The needles are sterile and individually packed.

COMFORTABLE INJECTION

European, premium needle quality, triple beveled tip and precision lancet point grinding.

THIN WALL TECHNOLOGY

Ensures optimal insulin flow.

SILICONIZED NEEDLE

Minimizes irritation and gives a smooth penetration.

SHORT NEEDLES

Reduce risk of a painful injection into the muscle.

VERSATILE

Fit most commonly used insulin pens.

COLOUR-CODING

Simple identification of different needle lengths.

FDA CLEARED

Technical Specifications

Colour Code		Specification	Article No.
Yellow		29G x 10 mm	CHPN2910
Pink		29G x 12.7 mm	CHPN29127
Dark blue		30G x 6 mm	CHPN3006
Light blue		30G x 8 mm	CHPN3008
Green		31G x 4 mm	CHPN3104
Dark violet		31G x 5 mm	CHPN3105
Dark blue		31G x 6 mm	CHPN3106
Light blue		31G x 8 mm	CHPN3108
Green		32G x 4 mm	CHPN3204
Dark violet		32G x 5 mm	CHPN3205
Dark blue		32G x 6 mm	CHPN3206
Green		33G x 4 mm	CHPN3304 *

MEDOJECT Fine Pen Needles are available in many sizes.
For more information Contact Sales Department.
* Not FDA cleared.

Shelf life 5 years

Packaging:
Primary packaging: 100 pcs.
Transit packaging: 2,000 pcs.
Pallet: 160,000 pcs.





Vitrex Safety Pen Needles With Dual Automatic Protective Shields

Vitrex Safety Pen Needles are intended for a safe and simple insulin injection via most insulin pens. The application is the same as normal insulin pens, but they offer an automatic double protection without extra operations after using. The coloured indicator shields in the front and back-end provide a visual, physical and audible confirmation of safety. The needles are sterile and individually packed.

SAFETY FUNCTION

Coloured automatic safety shields on front and back-end providing double protection after using.

COMFORTABLE INJECTION

Triple sharpened needle and top lubrication technology, which offers minimum penetration force.

SIMPLE AND SAFE APPLICATION

Easy to trigger the safety function after injection. Same application as normal pen needles, but with a safety function both blocking top and bottom needle after use.

VERSATILE


Uniform hub screw which fits most commonly used insulin pens.

Technical Specifications


Colour Code		Specification	Article No.
Purple		30G x 5 mm	531113
Dark blue		30G x 8 mm	531213
Yellow		31G x 5 mm	531313
Orange		31G x 6 mm	531413
Green		32G x 4 mm	531513

Shelf life 5 years

Before use:



After use:



Packaging:
Primary packaging: 100 pcs.
Transit packaging: 800 pcs.
Pallet: 57,600 pcs.



Chirana® Insulin and Tuberculin Syringes

Chirana Insulin and Tuberculin syringes are with integrated needle for single use. The needle cap is made of LATEX-FREE rubber gasket. The syringes are individually packed in boxes of 100 pcs. They are pyrogen-free and non-toxic.

USER CONVENIENCE: Needle permanently fixed to syringe body.

USER SAFETY: Integrated needle cap.

COMFORT: Grinded, precision lancet point and silicone treated tip for smooth penetration.

PATIENT SAFETY: Single-use only. Sterilised by ETO. Dead space is minimized to reduce waste and transmission of infectious agents. The dead space does not exceed the limits (EN ISO 8537:2016, art. 5.11.1).

FDA CLEARED *



Packaging:
Primary packaging: 100 pcs.
Secondary packaging: 2,000 pcs.
Pallet: 48,000 pcs.

* FDA clearance for Insulin syringes.

Technical Specifications

Type	Vol. ml.	Needle	Size Metric	Article No:
Tuberculin	1	29G x 1/2"	12.7 mm	CHTUB0129
Tuberculin	1	Luer		CHTUB01
Insulin U 40	1	29G x 1/2"	12.7 mm	CHINS4129
Insulin U 100	1	27G x 1/2"	12.7 mm	CHINS0127
Insulin U 100	1	29G x 5/16"	8 mm	CHINS01298
Insulin U 100	1	29G x 1/2"	12.7 mm	CHINS0129
Insulin U 100	1	30G x 1/2"	12.7 mm	CHINS0130
Insulin U 100	1	30G x 5/16"	8 mm	CHINS01308
Insulin U 100	1	31G x 4/16"	6 mm	CHINS01316
Insulin U 100	1	31G x 5/16"	8 mm	CHINS01318
Insulin U 100	1	Luer		CHINS01
Insulin U 100	0.5	29G x 1/2"	12.7 mm	CHINS00529
Insulin U 100	0.5	30G x 1/2"	12.7 mm	CHINS00530
Insulin U 100	0.5	30G x 5/16"	8 mm	CHINS005308
Insulin U 100	0.5	31G x 5/16"	8 mm	CHINS005318
Insulin U 100	0.5	31G x 4/16"	6 mm	CHINS005316
Insulin U 100	0.3	29G x 1/2"	12.7 mm	CHINS00329
Insulin U 100	0.3	30G x 1/2"	12.7 mm	CHINS00330
Insulin U 100	0.3	31G x 5/16"	8 mm	CHINS003318
Insulin U 100	0.3	31G x 4/16"	6 mm	CHINS003316

CHIRANA Syringes are available in many sizes.
For more information Contact Sales Department.

Shelf life 5 years

 CHIRANA injecta
Member of Vitrex Medical Group



Chirana® 2-part disposable syringes and Chirana® 2-part disposable syringes with needles

Chirana 2-part disposable syringes without and with needle are made of polypropylene/polyethylene. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. The syringes are available with Luer tip. The syringe with needle is available in other needle combinations.

USER CONVENIENCE: Ergonomic design, smooth flow, highly transparent barrel, tight, black extended graduation, easy to read, double stop lines and guaranteed 0 ml position.

USER SAFETY: Individually blister packed and the needle is side-packed in a single blister package.

COMFORT: Designed tip for best connection, double stop lines and scratch-free for glove manipulation.

PATIENT SAFETY: Single-use only. Sterilised by ETO.

Chirana 2-part disposable syringes

Vol. (ml)	Pcs./Box	Pcs./Carton	Article No.
2	100	2800	CH002L
5	100	1900	CH005L
10	100	1200	CH010L
20	80	960	CH020L

Shelf life 5 years

Chirana 2-part disposable syringes with needles

Specification	Pcs./Box	Pcs./Carton	Article No.
2 ml with needle 0.60 x 30 mm (23G x 1 1/4")	100	2800	CH002L0630
5 ml with needle 0.70 x 35 mm (22G x 1 2/5")	100	1900	CH005L0735
10 ml with needle 0.80 x 40 mm (21G x 1 1/2")	100	1200	CH010L0840
20 ml with needle 0.80 x 40 mm (21G x 1 1/2")	80	960	CH020L0840

Shelf life 5 years



Chirana® 3-part disposable syringes and Chirana® 3-part disposable syringes with needles

Chirana 3-part disposable syringes without and with needle are made of polypropylene, polyethylene and lubricated synthetic rubber. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. The syringes are available with Luer and Luer lock. The syringe with needle is available in other needle combinations.

USER CONVENIENCE: Ergonomic design, smooth flow, highly transparent barrel, tight, black extended graduation and easy to read.

USER SAFETY: Individually blister packed and the needle is side-packed in a single blister package.

COMFORT: Designed tip for best connection, scratch-free for glove manipulation and centrally positioned Luer lock tip.

PATIENT SAFETY: Single-use only. Sterilised by ETO. Low dead space volume for 1 ml.

Chirana 3-part disposable syringes

Vol. (ml)	Pcs./Box	Pcs./Carton	Article No. LUER	Article No. LUER LOCK
1	100	2000	CH03001L	CH03001LL
2	100	3000	CH03002L	CH03002LL
3	100	3000	CH03003L	CH03003LL
5	100	2400	CH03005L	CH03005LL
10	100	1600	CH03010L	CH03010LL
20	50	800	CH03020L	CH03020LL
30	50	400	CH03030L	CH03030LL
50		100	CH03050L	CH03050LL

Shelf life 5 years

Chirana 3-part disposable syringes with needles

Specification	Pcs./Box	Pcs./Carton	Article No.
2 ml with needle 0.60 x 30 mm (23G x 1 1/4")	100	3000	CH03002L0630
5 ml with needle 0.70 x 30 mm (22G x 1 1/4")	100	2400	CH03005L0730
10 ml with needle 0.80 x 40 mm (21G x 1 1/2")	100	1600	CH03010L0840
20 ml with needle 0.80 x 40 mm (21G x 1 1/2")	50	800	CH03020L0840
50 ml Luer with needle 1.20 x 40 mm (18G x 1 1/2")	50	450	CH03050L1240

Shelf life 5 years



Chirana® Perfusion 3-part syringes

Chirana Perfusion 3-part syringes are made of polypropylene and lubricated synthetic rubber. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. They are intended for perfusion applications and available in transparent, opaque and black and with aspiration needle and filter needle.

USER CONVENIENCE: Black graduation, easy to read and double graduation on barrel and piston. Compatible with perfusion pumps and dialysis units of world leader manufacturers.

PATIENT SAFETY: Single-use only. Sterilised by ETO.

Specification	Pcs./Carton	Article No.
20 ml transparent for dialysis machine use	100	CH03020PT
30 ml transparent for dialysis machine use	100	CH03030PT
50 ml transparent	100	CH03050PT
50 ml opaque	100	CH03050PO
50 ml transparent with needle 2 x 30 mm	100	CH030502030PT
50 ml transparent BB	100	CH03050PTB
50 ml opaque BB	100	CH03050POB
50 ml transparent BB with needle 2 x 30 mm	100	CH030502030PTB
50 ml opaque BB with needle 2 x 30 mm	100	CH030502030POB
50 ml black BB with needle 2 x 30 mm	100	CH030502030PBB
50 ml transparent F	100	CH03050PTF
50 ml opaque F	100	CH03050POF
50 ml transparent F with needle 1,80 x 40 mm	100	CH030501840PTF

Shelf life 5 years



Chirana® 3-part disposable Catheter syringes

Chirana 3-part disposable Catheter syringes are made of polypropylene and lubricated synthetic rubber. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. The syringes are intended for irrigation application and available in Catheter tip and with Luer adapter for wider application range.

USER CONVENIENCE: Ergonomic design, smooth flow, highly transparent barrel, tight, black extended graduation and easy to read.

COMFORT: Designed tip for best connection and scratch-free for glove manipulation.

PATIENT SAFETY: Single-use only. Sterilised by ETO.

Specification	Pcs./Box	Pcs./Carton	Article No.
50 ml Catheter		100	CH03050C
50 ml Catheter with Luer adapter		100	CH03050CLA
100 ml Catheter with Luer adapter	25	100	CH03100C
150 ml Catheter Janette®		50	CH03150CS

Shelf life 5 years



Medoject® Hypodermic needles

Colour Code		Size Metric	Size G/"	Article No.
Light yellow	●	0.30 x 12 mm	30G x 1/2"	CH30012
Grey	●	0.40 x 12 mm	27G x 1/2"	CH27012
Grey	●	0.40 x 20 mm	27G x 4/5"	CH27045
Grey	●	0.40 x 20 mm	27G x 3/4"	CH27034
Grey	●	0.40 x 40 mm	27G x 1 1/2"	CH27112
Brown	●	0.45 x 12 mm	26G x 1/2"	CH26012
Brown	●	0.45 x 16 mm	26G x 5/8"	CH26058
Brown	●	0.45 x 25 mm	26G x 1"	CH26100
Orange	●	0.50 x 16 mm	25G x 5/8"	CH25058
Orange	●	0.50 x 20 mm	25G x 4/5"	CH25045
Orange	●	0.50 x 25 mm	25G x 1"	CH25100
Orange	●	0.50 x 40 mm	25G x 1 1/2"	CH25112
Purple	●	0.55 x 25 mm	24G x 1"	CH24100
Purple	●	0.55 x 38 mm	24G x 1 1/2"	CH24112
Blue	●	0.60 x 16 mm	23G x 5/8"	CH23058
Blue	●	0.60 x 25 mm	23G x 1"	CH23100
Blue	●	0.60 x 30 mm	23G x 1 1/4"	CH23114
Blue	●	0.60 x 40 mm	23G x 1 1/2"	CH23112
Blue	●	0.60 x 60 mm	23G x 2 1/4"	CH23214
Black	●	0.70 x 25 mm	22G x 1"	CH22100
Black	●	0.70 x 30 mm	22G x 1 1/4"	CH22114
Black	●	0.70 x 35 mm	22G x 1 2/5"	CH22125
Black	●	0.70 x 40 mm	22G x 1 1/2"	CH22112
Black	●	0.70 x 50 mm	22G x 2"	CH22200
Green	●	0.80 x 16 mm	21G x 5/8"	CH21058
Green	●	0.80 x 25 mm	21G x 1"	CH21100
Green	●	0.80 x 30 mm	21G x 1 1/4"	CH21114
Green	●	0.80 x 40 mm	21G x 1 1/2"	CH21112
Green	●	0.80 x 50 mm	21G x 2"	CH21200
Yellow	●	0.90 x 25 mm	20G x 1"	CH20100
Yellow	●	0.90 x 30 mm	20G x 1 1/4"	CH20114
Yellow	●	0.90 x 40 mm	20G x 1 1/2"	CH20112
Yellow	●	0.90 x 50 mm	20G x 2"	CH20200
Ivory	○	1.10 x 25 mm	19G x 1"	CH19100
Ivory	○	1.10 x 30 mm	19G x 1 1/4"	CH19114
Ivory	○	1.10 x 40 mm	19G x 1 1/2"	CH19112
Ivory	○	1.10 x 50 mm	19G x 2"	CH19200
Pink	●	1.20 x 25 mm	18G x 1"	CH18100
Pink	●	1.20 x 40 mm	18G x 1 1/2"	CH18112
Pink	●	1.20 x 50 mm	18G x 2"	CH18200
Red purple	●	1.50 x 50 mm	17G x 2"	CH17200
White	○	1.60 x 25 mm	16G x 1"	CH16100
White	○	1.60 x 40 mm	16G x 1 1/2"	CH16112
Blue grey	●	1.80 x 40 mm	15G x 1 1/2"	CH15112

Shelf life 5 years

Medoject®

Medoject needles are made of polypropylene hub and stainless steel tube. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. The needles are intended for intramuscular, subcutaneous, intravenous and intra-arterial application. They are individually packed in boxes of 100 pcs, cartons of 4000 pcs.

USER CONVENIENCE: Colour-coded transparent hub for flow visual control, metric and gauge measurements on packaging and mainly thin wall tube material.

USER SAFETY: Safety needle is available.

COMFORT: Grinded, precision lancet point and silicone treated for smooth penetration.

PATIENT SAFETY: Single-use only. Sterilised by ETO.

Medoject® Hypodermic needles

Medoject Hypodermic needles are intended for intramuscular, subcutaneous, intravenous and intra-arterial application. There is possibility of short bevel lancet point grinding.



Medoject® Safety Hypodermic needles

Medoject Safety Hypodermic needles have an integrated safety feature for one-handed manipulation.

Colour Code		Size Metric	Size G/"	Article No.
Blue	●	0.60 x 30 mm	23G x 1 1/4"	CHS23114
Black	●	0.70 x 40 mm	22G x 1 1/2"	CHS22112
Green	●	0.80 x 40 mm	21G x 1 1/2"	CHS21112
Yellow	●	0.90 x 40 mm	20G x 1 1/2"	CHS20112

Shelf life 5 years

Medoject® Blunt fill needles

Medoject Blunt fill needles are also available with filter of 5µm for aspiration filtration procedure, which risks particle elimination during procedure.

Colour Code		Size Metric	Size G/"	Article No.
Red	●	1.20 x 40 mm Blunt fill	18G x 1 1/2"	CH18112BF
Pink	●	1.20 x 40 mm Blunt fill w/ filter	18G x 1 1/2"	CH18112F
Pink	●	1.20 x 50 mm Blunt fill w/ filter	18G x 2"	CH18200F

Shelf life 5 years



Medoject® Special Hypodermic needles

Medoject Special Hypodermic needles are intended for a wide special application range. There is possibility of short bevel lancet point grinding.

Colour Code		Size Metric	Size G/"	Article No.
Blue	●	0.60 x 80 mm	23G x 3 1/8"	CH23318
Green	●	0.80 x 100 mm	21G x 4"	CH21400
Green	●	0.80 x 120 mm	21G x 4 3/4"	CH21434
Yellow	●	0.90 x 70 mm	20G x 2 4/5"	CH20245
Pink	●	1.20 x 70 mm	18G x 2 4/5"	CH18245
Pink	●	1.20 x 100 mm	18G x 4"	CH18400
White	○	1.60 x 70 mm	16G x 2 4/5"	CH16245
Blue grey	●	1.80 x 70 mm	15G x 2 4/5"	CH15245
Pale green	●	2.00 x 70 mm	14G x 2 4/5"	CH14245
Pale green	●	2.00 x 80 mm	14G x 3 1/8"	CH14318
Pale green	●	2.00 x 100 mm	14G x 4"	CH14400
Pale green	●	2.00 x 120 mm	14G x 4 3/4"	CH14434

Other sizes available upon request.

Shelf life 5 years

Inox® Ophthalmic needles

Inox Ophthalmic needles has a blunt end for ophthalmic application, which includes scratch-free edges for comfort use.

Colour Code		Size Metric	Size G/"	Article No. STRAIGHT	Article No. BENT
Grey	●	0.40 x 20 mm	27G x 4/5"	CHOP27045S	CHOP27045B
Blue	●	0.60 x 20 mm	23G x 4/5"	CHOP23045S	CHOP23045B

Shelf life 5 years



Chiraflex® IV Cannulas and Chiraflex® Safety IV Cannulas

Chiraflex IV Cannulas and Chiraflex Safety IV Cannulas are intended for peripheral vascular application. The catheter is made from Teflon material and available in Safety. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. They are packed in boxes of 50 pcs, cartons of 1000 pcs.

USER CONVENIENCE: Compatible with all standard medical devices for wide application.

USER SAFETY: Safety cannula is available.

PATIENT SAFETY: Single-use only. Sterilised by ETO.

Colour Code		Size G/''	Size Metric	Article No. With Port With Wings	Article No. Without Port With Wings	Article No. Without Port Without Wings	Article No. SAFETY With Port With Wings
Purple	●	26G	0.60 x 19	CHFX0102619	CHFX0202619	CHFX0302619	CHFXS0102619
Yellow	●	24G	0.70 x 19	CHFX0102419	CHFX0202419	CHFX0302419	CHFXS0102419
Blue	●	22G	0.90 x 25	CHFX0102225	CHFX0202225	CHFX0302225	CHFXS0102225
Pink	●	20G	1.10 x 32	CHFX0102032	CHFX0202032	CHFX0302032	CHFXS0102032
Green	●	18G	1.30 x 45	CHFX0101845	CHFX0201845	CHFX0301845	CHFXS0101845
White	○	17G	1.50 x 45	CHFX0101745	CHFX0201745	CHFX0301745	CHFXS0101745
Grey	●	16G	1.70 x 45	CHFX0101645	CHFX0201645	CHFX0301645	CHFXS0101645
Orange	●	14G	2.10 x 45	CHFX0101445	CHFX0201445	CHFX0301445	CHFXS0101445

Shelf life 5 years



Chiraflex® Mandrins

Chiraflex Mandrins are intended for longer interruptions in infusion and has a screw onto IV cannula. The mandrins are compatible with Chiraflex IV Cannulas and Chiraflex Safety IV Cannulas. They are colour-coded as per Chiraflex IV Cannula colour-coding.

Colour Code	Size	Pcs./Box	Pcs./Carton	Article No.
Blue	22G	50	1000	CHFX20022
Pink	20G	50	1000	CHFX20020
Green	18G	50	1000	CHFX20018
White	17G	50	1000	CHFX20017
Grey	16G	50	1000	CHFX20016
Orange	14G	50	1000	CHFX20014

Shelf life 5 years

Chiraflex® Scalp vein sets

Chiraflex Scalp vein sets are intended for shorter time infusion. The set includes a stainless steel needle with a length of 20 mm and a 30 cm plastic tube.

Colour Code	Size	Diameter (mm)	Pcs./Box	Pcs./Carton	Article No.
Grey	27G	0.40	100	2000	CH260001127
Orange	25G	0.50	100	2000	CH260001125
Blue	23G	0.60	100	2000	CH260001123
Black	22G	0.70	100	2000	CH260001122
Green	21G	0.80	100	2000	CH260001121
Yellow	20G	0.90	100	2000	CH260001120
Ivory	19G	1.10	100	2000	CH260001119
Pink	18G	1.20	100	2000	CH260001118

Shelf life 5 years





Chiraplus™ Infusion set

Chiraplus Infusion sets include 15 µm filter, ventilation port and Luer lock adapter. Chiraplus P is intended for pressure infusion and does not include flash ball. Chiraplus G is intended for gravity infusion and includes flash ball.

Product	Chamber Vol. (ml)	Tube length (mm)	Packaging Type	Pcs./Carton	Article No.
Chiraplus P	10.5	1500	Blister	100	CHP336
Chiraplus P Dehp free	10.5	1500	Blister	100	CHP536
Chiraplus P 12 Dehp free	12	1500	Blister	100	CHP506
Chiraplus P 190	10.5	1900	Blister	100	CHP339
Chiraplus P 190 Safety	10.5	1900	Blister	100	CHP338
Chiraplus P Dehp free Autofusion 180	15	1800	Blister	300	CHP576
Chiraplus G	10.5	1500	Blister	100	CHG331
Chiraplus G	10.5	1500	Poly bag	500	CHG317
Chiraplus G Dehp free	10.5	1500	Blister	100	CHG531
Chiraplus G 12	12	1500	Blister	100	CHG301
Chiraplus G 12 Dehp free	12	1500	Blister	100	CHG501
Chiraplus G	8.5	1500	Poly bag	500	CHG318
Chiraplus G Dehp free	8.5	1500	Poly bag	500	CHG518
Chiraplus G with needle	8.5	1500	Poly bag	500	CHG319
Chiraplus G with metal spike	8.5	1500	Poly bag	500	CHG316
Chiraplus G PVC free	12	1500	Blister	100	CHG731
Chiraplus G PVC free opaque	12	1900	Blister	100	CHG901

Shelf life 5 years



Chirahem™ Transfusion set

Chirahem Transfusion set is intended for blood transfusion administration. The set has a 200 µm blood filter.

Product	Chamber Vol. (ml)	Tube length (mm)	Packaging Type	Pcs./Carton	Article No.
Chirahem	17	1500	Blister	100	CHT314
Chirahem Dehp free	17	1500	Blister	100	CHT514
Chirahem	17	1500	Poly bag	300	CHT320
Chirahem with needle	17	1500	Poly bag	300	CHT321
Chirahem with metal spike	17	1500	Poly bag	300	CHT322
Chirahem latex free	17	1500	Blister	100	CHT518

Shelf life 5 years



Accessories IV Administration

Within IV Administration for Chiraplust Infusion set and Chirahem Transfusion set we offer following accessories:

Chiraplust™ Spikes: Intended for comfortable and safe drug manipulation.

Chiraplust™ Injection ports: Offer additional drug administration when applied with other devices.

Chiraplust™ Stoppers: Give a short-term option for storage or transfer of solution.

Chiraline™ Extension lines: Extension tubing in 45 cm or 180 cm, where Luer lock ends are available.

Chiraway™ 3-way stopcock: Available in standard or lipid resistant material with colour-coding differentiation.

Product	Pcs./Box	Pcs./Carton	Article No.
Infusion spike Red air 0.1 µm + 5 µm solution filter	100	1000	CHIS01
Infusion spike Blue air 0.1 µm + 5 µm solution filter	100	1000	CHIS02
Infusion spike Green air 0.1 µm filter	100	1000	CHIS03
Needle free injection port	50	2500	CHNF01
Injection port	100	5000	CHIP01
Combistopper Red	100	5000	CHCS01
Combistopper Blue	100	5000	CHCS02
Combistopper White	100	5000	CHCS03
Stopper Luer lock	100	5000	CHLL01
Chiraline 450	120	960	CHL318045
Chiraline 1800	85	680	CHL318180
Chiraway Standard	50	1000	CHW001
Chiraway Lipid resistant	50	1000	CHW002

Shelf life 5 years



Examination

Chiraskin™ Gloves

Chiraskin Gloves are non-sterile and intended for examination. They offer high elasticity, high strength, excellent grip sensitivity, safe protection and good touch sense.

Size	Pcs./Box	Pcs./Carton	Article No. Latex Powdered	Article No. Latex Powder Free	Article No. Nitrile Powder Free
Gloves Small	100	1000	CHL01S	CHLF01S	CHNF01S
Gloves Medium	100	1000	CHL01M	CHLF01M	CHNF01M
Gloves Large	100	1000	CHL01L	CHLF01L	CHNF01L

Shelf life 5 years

Examination products


The examination products are intended for variety of application range to support safe and clean healthcare.

Product	Pcs./Box	Pcs./Carton	Article No.
Face mask Tie On, 3-ply, non-woven	50	1000	CHFM01
Face mask Ear Loop, 3-ply, non-woven	50	1000	CHFM02
Face mask Ear Loop, 3 ply, Type IIR	50	2500	V81308
Bouffant cap 21", non-woven	100	1000	CHC01
Bouffant cap 24", non-woven	100	1000	CHC02
Shoe covers, CPE, blue	100	2000	CHSC01
Bed cover, CPE, blue	10	100	CHBC01

Shelf life 5 years



CHIRMAX surgical
Member of Vitrex Medical Group





When experience matters

Our surgical sutures are manufactured in line with state of the art technology. Each individual production process is subject to strict scientific quality control.

Continuous inspection of both mechanical and microbiological characteristics is given top priority. We have more than 90 years of experience in manufacturing surgical sutures.

Our team understands the specific needs of medical professionals by providing them with access to the most complete portfolio of surgical sutures. The result is the current modern product range of absorbable and non-absorbable sutures, which consists of more than 2500 individual surgical suture products.

Production takes place in accordance with GMP guidelines, and a comprehensive quality management system in line with ISO 13485 has been implemented.

All products have been CE marked since 1999.



Surgical Meshes

Surgical Meshes have become standard for the repair of abdominal wall hernias. Surgical mesh is a sterile woven material designed for permanent implantation within the body during open or laparoscopic procedures.

A wide range of mesh implants are available, with two main functions: to stabilize and strengthen soft tissue defects, and to act as a sling to support prolapsed organs and viscera. It is most commonly used in pelvic organ prolapse, urinary stress incontinence and hernia repairs.

CAPROMESH meshes

Partially absorbable

The sterile, partially absorbable, composite mesh is made from two components: one absorbable and one non-absorbable. After absorption of the absorbable part, only the polypropylene mesh remains in the body.

Indication: Capromesh is indicated in the operative treatment of weakened abdominal wall and fascial defects - in particular hernias, the healing of which requires long-term bridging or strengthening of the damaged structures.

CHIRALEN meshes

Non-absorbable

The sterile and surgical meshes are made from undyed polypropylene filaments. The knitted filaments have a smooth and hydrophobic surface, which biologically are absolutely inert. Chiralen meshes are available in different sizes, all offering high strength and high flexibility. Chiralen meshes are available in a weight of approx. 110 g/m² or in the light version with a weight of approx. 30 g/m².


Indication: Chiralen meshes are indicated in the operative treatment of fascial defects, in particular hernias, which require long-term bridging of damaged structures.



Chiralen Mesh



Capromesh



Chiralen Mesh Light

Surgical Sutures - Absorbable threads

CHIRASORB rapid braided

Intended for use in Skin closure, Paediatric surgery, Episiotomies, Circumcision, Closure of oral mucosa and Ophthalmic surgery for conjunctival sutures.



Composition	Colour	Range	Characteristics
Braided and coated suture made from 90% glycolide and 10% L-lactide	Violet and Natural	6/0 to 2 USP	Absorbable, multifilament, coated, synthetic, high tensile strength, reduced capillary effect, short-term healing

CHIRLAC rapid braided

Intended for use in Gynaecology, Skin Closure, Ophthalmic, Odontology and Paediatric surgery.



Composition	Colour	Range	Characteristics
Braided and coated suture made of Polyglycolic acid (PGA)	Violet and Natural	6/0 to 4 USP	Absorbable, multifilament, coated, synthetic, high tensile strength, reduced capillary effect, short-term healing

MONOLAC monofilament

General soft tissue approximation and/or ligation in all specialties, where an absorbable material is indicated.



Composition	Colour	Range	Characteristics
Monofilament suture made of Glycolide and ϵ -caprolactone	Violet and Natural	6/0 to 2 USP	Absorbable, monofilament, synthetic, high tensile strength, excellent handling properties, mid-term healing

CHIRASORB braided

Intended for use in general soft tissue approximation and/or ligation in all specialties, including use in Ophthalmic surgery.



Composition	Colour	Range	Characteristics
Braided and coated suture made from 90% glycolide and 10% L-lactide	Violet and Natural	10/0 to 6 USP	Absorbable, multifilament, coated, synthetic, high tensile strength, reduced capillary effect, mid-term healing

CHIRASORB Plus braided

Antibacterial

Intended for use in general soft tissue approximation and/or ligation in all specialties.



Composition	Colour	Range	Characteristics
Antibacterial, braided and coated suture made from 90% glycolide, 10% L-lactide and chlorhexidine diacetate	Violet and Natural	5/0 to 4 USP	Absorbable, multifilament, anti-bacterial coating, synthetic, high tensile strength, reduced capillary effect, mid-term healing

CHIRLAC braided


Intended for use in general soft tissue approximation and/or ligation, surgical adaptation of Muscles and Fascia, suturing of Parenchymal organs, Peripheral nerves and ligatures of small vessels and for Intradermal suturing.




Composition	Colour	Range	Characteristics
Braided and coated suture made of Polyglycolic acid (PGA)	Violet and Natural	9/0 to 6 USP	Absorbable, multifilament, coated, synthetic, high tensile strength, reduced capillary effect, mid-term healing

Polydiox monofilament

Intended for use in general soft tissue approximation, Orthopaedics, Fascia Closure, Urology, Gynaecology and Plastic surgery. Recommended for long-term healing.




Composition	Colour	Range	Characteristics
Monofilament suture made of Poly-p-dioxanone	Violet and Natural	7/0 to 4 USP	Absorbable, monofilament, synthetic, high tensile strength, excellent knotting properties, long-term healing



Surgical Sutures - Non-absorbable threads

CHIRAFLON monofilament

Intended for use in all surgical specialities, especially Orthopaedic, Traumatology (tendons), Plastic and Ophthalmic surgery (optimal, temporary, easily, removable suture).



Composition	Colour	Range	Characteristics
Monofilament suture made of Polyvinylidene Fluoride (PVDF) - similar to Teflon	Blue	8/0 to 3 USP	Non-absorbable, monofilament, synthetic, no capillary effect, very smooth, excellent knotting properties, easy to handle, physiologically inert

Chiralen monofilament


Intended for use in all surgical specialities, especially Cardiovascular, Plastic, Reconstructive, Ophthalmic, Orthopaedics and Traumatology (tendons).



Composition	Colour	Range	Characteristics
Monofilament suture made of Polypropylene	Blue and Black	10/0 to 2 USP	Non-absorbable, monofilament, synthetic, no capillary effect, smooth, non-porous surface, suitable for stitches in infected wounds

SILON monofilament

Intended for use in Plastic and General surgery and wherever suturing the skin.



Composition	Colour	Range	Characteristics
Monofilament suture made of Polyamide 6	Blue and Black	10/0 to 2 USP	Non-absorbable, monofilament, synthetic, very smooth, excellent knotting properties, easy to handle, loss of tensile strength after long-term implantation

TERVALON braided

Intended for use whenever non-absorbability and low friction passage through tissues are needed. It is especially useful in suturing Tendons, Schwann sheaths and skin in Plastic and Reconstructive surgery.



Composition	Colour	Range	Characteristics
Braided and coated suture made of Polyester	Green and White	6/0 to 6 USP	Non-absorbable, multifilament, synthetic, high tensile strength, no capillary effect, well tolerated by the tissue, low friction

SILON braided

Intended for use in Skin Closure, Subcutaneous tissue, Transversely cut muscles, Fascia, Peritoneum, Gastrointestinal Tract and Blood vessels ligation.



Composition	Colour	Range	Characteristics
Braided suture made of Polyamide	White and Black	5/0 to 5 USP	Non-absorbable, multifilament, synthetic, high tensile strength, excellent knotting properties, loss of tensile strength after long-term implantation


Silk braided

Intended for use in Skin closure, Subcutaneous tissue, Fascia, Peritoneum, Plastic surgery, Gynaecology and the organs of the Gastrointestinal tract.



Composition	Colour	Range	Characteristics
Braided and coated suture made of Silk	Black	9/0 to 5 USP	Non-absorbable, multifilament, synthetic, no capillary effect, loss of tensile strength after long-term implantation

For more information about sutures and meshes, and catalogue of full product range, please contact us.







CHIRANA T. Injecta, a.s.
Nám. Dr. Alberta Schweitzera 194
916 01 Stará Turá, Slovakia
Phone: +421 32 7709 901

Vitrex Medical A/S
Vasekaer 6-8
2730 Herlev, Denmark
Phone: +45 4494 7011

CHIRANA T. Injecta, s.r.o.
Komořanská 2148
143 00 Praha 4
Czech Republic
Phone: +420 241 097 904



Vitrex Medical Group 2024-11 Rev.001

